

Portable Scale



- Ideal for laboratory, industrial, education and home use; lightweight
- **Applications:** forensic testing, quality control, formulation, soil sampling, postal weighing & dietary needs
- **Units of measure:** g, lb:oz, Newton
- Superior RF protection for use near cell phones and production areas
- **Features:** mechanical/software overload protection, low battery indicator, auto shut-off, push button calibration, large LCD
- **Includes:** large, removable stainless steel pan; AC Adapter or batteries
- ABS Compact housing, low profile

| | M102299 | M102300 | M102301 |
|----------------|-------------------------------------|--------------------------------------|--------------------------------------|
| Capacity: | 200 g | 0.1 g | 145 x 133 mm |
| Readability: | 2000 g | 1 g | 145 x 133 mm |
| Platform Dia.: | 5000 g | 1 g | 145 x 133 mm |
| Cal. Masses: | 200 g Calibration Mass (M102317) | 1000 g Calibration Mass (M102318) | 2000 g Calibration Mass (M102319) |
| M288402: | Weighing Bowl, 1050 ml | | |



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002