


Scale, Portable, Navigator™



- Ideal for a QC, inventory, warehouse, food, portioning or lab
- Superior mechanical overload protection system (bears loads up to 4X its rated capacity) & ultra fast response time (1 sec.)
- Touchless; wave of hand controls tare, zero, display & print also avoids cross-contamination
- Two sensors can be set up independently to control up to 36 combinations of operations
- **Applications:** counting, % weighing, checkweighing, display hold and accumulation/totalization
- Stainless steel pan (9.45 x 9.82 in.) & a corrosion-free plastic housing has a spill-resistant design that is easy to clean
- **Units of measure:** g, kg, lb, oz, lbs.:oz, troy oz, grain, Newton Carat, Pennywght, teal, ticl, tola
- **Includes:** internal rechargeable battery (100 hrs); AC adapter
- USB, RS232 or Ethernet Option
- **Option:** full-body in-use cover, carry case • 2 YR WARRANTY

	Cap. x Read.	Repeatability	Linearity
M110965	1,600 g x 0.1 g	0.1 g	± 0.2 g
M110966	3,200 g x 0.2 g	0.2 g	± 0.4 g
M110967	6,400 g x 0.5 g	0.5 g	± 1 g
M110968	10,000 g x 0.5 g	0.5 g	± 1 g
M110969	16,000 g x 1 g	1 g	± 2 g



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002