

Personal Water Bath, 4L MyBath™



- Compact, 4L bath; less than 8" wide
- Digital accuracy and convenience
- Quik-CAL™ user calibration function
- Hinged lid provides a secure, covered environment; can be flipped open (to rear) or completely removed
- **Includes:** tube racks for 1.5/2.0ml, 15ml and 50ml tubes
- 2 YEAR WARRANTY



	M111455
Capacity:	4L
Range:	Ambient 5 to 100°C
Accuracy:	±0.2°C
Increments:	0.1°C
Uniformity:	±0.2°C
Chamber (WxDxH):	5.3 x 9.4 x 5.5 in. / Stainless Steel
Overall Dims. (WxDxH):	7.7 x 13.9 x 9 in. / H: 12.2 in. w/ lid
Weight:	15 lbs.

Cat. No.	Description
M111455	4L MyBath™
M111456	Optional Test Tube Rack for 40 x .05ml Tubes