

# Large Capacity Incubators



- Compact, low-cost and dependable
- Designed for general purpose as well as clinical, industrial, school laboratories or physician's offices
- Has a full PID microprocessor control that maintains accurate temperature
- Settings within 1°C, even in varying ambient conditions
- Has a unique temperature tracking feature that stores maximum process temperature deviation from set point in 0.1°C units
- Cold rolled steel-ivory baked-on enamel exterior
- Interior is aluminum and fiber glass insulated with chrome plated adjustable shelf and fixed bottom shelf
- See-thru tinted acrylic door for viewing all samples
- Beakers, flasks, petri dishes and test tube trays can be processed together
- A lab thermometer is included
- Additional wire shelves with adjustable supports are available
- Each incubator is tested and approved before shipment
- °C or °F readout
- 2 YEAR WARRANTY (on heating element)

|                         | M279505          | M279502          | M107590          | M107591          |
|-------------------------|------------------|------------------|------------------|------------------|
| Stability at 35°C:      | ± 1/2 °C         |                  |                  |                  |
| Display Resolution:     | 1.0° or 0.1°     |                  |                  |                  |
| Temperature Uniformity: | Ambient to 62°C  |                  |                  |                  |
| Dimensions:             | 12 x 10 x 10 in. | 18 x 16 x 12 in. | 12 x 10 x 10 in. | 18 x 16 x 12 in. |



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002