

Benchtop pH Meter

Orion Star™



- Large, informative LCD screen clearly shows information, including pH or mV reading, temperature, meter mode and calibration
- AUTO-READ; locks in the stable reading on screen and ready indicator alerts when readings are stable
- Up to 3 point pH calibration
- Automatic recognition of USA/NIST and DIN buffers
- Non-volatile memory holds up to 50 data points
- Automatic recognition of USA/NIST and DIN buffers
- Easily recall calibration slope data for procedure checks
- Electrode arm (included) and probe holder make it easier to maintain and place probes into samples
- Power adapter (included), or 4 AA batteries (sold separately)
- The IP54-rated housing can handle splashes and can be wall-mounted if table space is a problem
- 3 YEAR WARRANTY ON METER

M111356 Meter Kit includes: 3-in-1 pH/ATC probe; Solution kit (60mL each of pH 4, 7, 10 buffers; storage solution & cleaning solution); electrode arm with redesigned holder and power adapter.

M111355		
pH	Range:	-2.000 to 16.00
	Resol.:	0.1, 0.01
	Relative Accuracy:	±0.01
Temp.	Range:	23° to 221°F (-5° to 105°C)
	Resol.:	0.1
	Relative Accuracy:	±0.1°C
mV	Range:	mV: ±1600.0; RmV: ±199.9 mV
	Resol.:	0.1
	Relative Accuracy:	±0.2 mV or 0.05% (greater)



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002