

## Balance, Portable Scoutn Pro®



- Ideal for laboratory, industrial or education use
- Superior shock resistance & overload protection protects wtgh cell
- **Applications:** Parts Counting, % Weighing, Totalization & Display Hold
- **Units of measure:** g, kg, oz, lbs., lbs.:oz, oz t, dwt
- **Features:** Removable platform, high contrast LCD, security bracket, shipping lock, weigh below hook, stability indicator, over/underload indicators, low battery indicator, auto shut-off
- Sealed front panel & molded spill ring to protect from spills
- Calibration mass included on 200g, 400g and 600g models
- **Includes:** stainless steel pan; AC Adapter or batteries
- Option of either RS232 or USB connectivity
- **Optional Accessories:** RS232 Interface, Hard Shell Carrying Case, Sample Scoops and Calibration Masses
- 2 YEAR WARRANTY

		Capacity (g)	Readability (g)	Repeatability (g)	Linear (g)
<b>Pan Size: 4.7 in. Dia.:</b>					
	<b>M102938:</b>	200	0.01	0.01	±0.01
	<b>M106589:</b>	400	0.01	0.01	±0.01
	<b>M102939:</b>	400	0.1	0.1	±0.1
<b>Pan Size: 6.5 x 5.5 in. (WxD):</b>					
	<b>M102940:</b>	600	0.1	0.1	±0.1
	<b>M102943:</b>	2000	0.1	0.1	±0.1
	<b>M106591:</b>	4000	0.1	0.1	±0.1
	<b>M102944:</b>	6000	1.0	1.0	±1.0
<b>Accessories:</b>	<b>M288402:</b>	Security Device			
	<b>M107745:</b>	RS232 Interface			
	<b>M288400:</b>	Hard Shell Case			
	<b>M107746:</b>	USB Interface			



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002