

Commercial / Industrial Balance



- Ideal for industrial and commercial food service; kitchens, salad bars, portion control, testing warehouse & production areas, QC & hydroponics
- Electronic, large LCD, easy to navigate
- **Indicators:** low battery, scale at zero, reading is stable
- **Units of measure:** g, kg, oz, lbs., lbs.:oz
- **Features:** All Stainless Steel body and removable platter, counting feature (secondary), Zero/Tare key, selectable auto backlight, auto off
- Adapter (incl.) or batteries (not incl.)
- 1 YEAR WARRANTY

| | M110954 |
|--------------------|--------------|
| Capacity: | 5000 g |
| Graduation: | 1 g |
| Platform Diameter: | 170 x 180 mm |



www.LABMART.com • 800-684-1234 • 908-561-1234 • Fax: 908-561-3002